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Nitrile-Examination Gloves – Instruction for use

Gloves for single use compliant with PPE Regulation (EU) 2016/425, EN ISO 21420:2020, EN ISO 374-1:2016 + A1:2018 and EN ISO 374-5:2016

Instructions for use are to be used in combination with the specific product-related information on each product packaging. The gloves are sold as a bundled unit with the shipping carton. In case products in the bundle are unbundled and sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit. The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonized European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms, and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.

EXPLANATION OF THE STANDARDS AND PICTOGRAMS							
EN ISO 374-1:2016 + A1:2018	Performance levels against permeation are defined based on the breakthrough times (tested according to EN 16523-1:2015+A1:2018) as follows:						
Type A / B / C	Performance level against permeation according to EN ISO 374-1:2016 + A1:2018	1	2	3	4	5	6
	Breakthrough time in minutes	>10	>30	>60	>120	>240	>480
 ABCDEFGHIJKLMNPST	Type A = Breakthrough time of chemical > 30 minutes against at least 6 chemicals on the list. Type B = Breakthrough time of chemical > 30 minutes against at least 3 chemicals on the list. Type C = Breakthrough time of chemical > 10 minutes against at least 1 chemical on the list. TEST CHEMICALS: A = methanol / B = acetone / C = acetonitrile / D = dichloromethane / E = carbon disulfide / F = toluene / G = diethylamine / H = tetrahydrofuran / I = ethyl acetate / J = n-heptane / K = sodium hydroxide 40 % / L = sulfuric acid 96 % / M = nitric acid 65 % / N = acetic acid 99 % / O = ammonia water 25 % / P = hydrogen peroxide 30 % / S = hydrofluoric acid 40 % / T = formaldehyde 37 %						
EN ISO 374-4:2019	The degradation (in %) indicates the change in the puncture resistance of the gloves after exposure to the respective test chemical.						

This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions using samples taken from the palm only and relates only to the chemical gloves suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemicals in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation resistance can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfection.

EN ISO 374-5: 2016  Virus	Tested for permeability according to EN ISO 374-2:2019 Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5:2016 / ASTM F1671 (virus resistance). Resistance to bacteria and fungi – pass Resistance to viruses – pass The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.						
 XXXX	XXXX = Notified Body responsible for EU type examination and supervising ongoing conformity						
	Gloves for single-use only and non-reusable.		Read the instructions for use carefully before use.		Temperature delimitation		
	Protect from moisture		Protect from sunlight		This is a non-sterile product.		

Precautions for use
 Before use, always check the gloves for possible mechanical damage, e.g. holes or tears. Do not use damaged gloves. Select gloves with a length appropriate for you, as the risk to the wrist area is minimal. Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers or others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. When removing the gloves, avoid the contact of gloves' outer surface with your skin, because the surface may have been contaminated with blood or other body fluids. Avoid snagging, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.

Ingredients/Hazardous Components
 Certain gloves may contain ingredients known to cause skin irritation or allergic reactions in sensitive individuals. Check the warnings on the specific packaging. Formulation will be shared upon request.

Storage Instructions
 The storage area should be cool, dry, and dust-free. Protect gloves from UV light sources, sunlight, oxidizing agents, and ozone. Store in the original packaging according to the temperature limits indicated on the packaging.

Disposal Instructions
 Dispose of gloves according to the valid regulations for the glove material. Gloves contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.

Alto Products Im- und Export GmbH, Paul-Thomas-Straße 48, 40599 Dusseldorf, Germany, ww.ahs-alto.de 01/2024

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